

With Technology Relieve Fatigue

Long time manual operation of root canal preparation will bring hand fatigue, especially the continuous preparation of multiple root canals or preparation of root canals with complex shapes. Endo motor can increase the efficiency a lot, to shorten time of preparation, so as to enable dentists' easier root canal preparations.



The choice of millions of dentists worldwide


Scan and Login website
for more information



Guilin Woodpecker Medical Instrument Co., Ltd.
Information Industrial Park, Guilin National High-Tech
Zone, Guilin, Guangxi, 541004 P. R. China

Tel:
Europe Sales Dept.: +86-773-5873196, +86-773-2125222
North America, South America &
Oceania Sales Dept.: +86-773-5873198, +86-773-2125123
Asia & Africa Sales Dept.: +86-773-5855350, +86-773-2125896
Fax: +86-773-5822450
E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com
Website: <http://www.glwoodpecker.com>
ZMN/WI-10-073 V1.0

DISTRIBUTED BY,



ENDO RADAR

Endo Motor



Intelligent touching control

Intelligent touching pen is easy to clean
and endows comfortable hand feeling.
Luminescent screen helps
easily identifying keys.



LCD Screen

Clear display and information feedback.



Continuous mode

Three automatic reverse functions,
avoiding instrument fractures result from
improper force exertion in manual operation.
When reaching the preset torque,
automatically reversely rotate; it will back to
forward rotation after disappearance of resistance.
When reaching the preset torque,
automatically reversely rotate; it will
automatically stop after disappearance of resistance.
Automatically stop when reaching preset torque.
Torque: 0.6-4.0Ncm
Speed: 100-650RPM



Reciprocating mode

-Be compatible with all popular
single file systems on the market.

Accurate angle control;
three optional reciprocating angles.



Program mode

Eight optional self-define modes
which can be automatically stored after setting.



Simplification

Control handpiece with wireless signal,
eliminating constraint of cord
and enabling flexible operation.
Light handpiece brings better handling feeling.



Long-endurance

Optional wireless charging and wire charging
Built-in battery in base realizes charging
by easily putting handpiece into the base.
Non-contact charging without external contact
point makes it not easy to age.

