

EASE-01D DENTAL UNIT

Led light, hand control and induction dual control system provide more convenient operation

Both strong and weak suction have a delay function to avoid pipeline blockage

The micro switches are moved forward to the front, it is more convenient to operate

The hanger can be rotated 90 degrees to adjust the position according to the doctor's habits

Revision: 202302
Youjoy reserves the right to make modifications without prior notice.



Ningbo Youjoy Dental Equipment Co.,Ltd.
Add: No.275,Lane321,Tongning Road,Jiangbei Investment & Pioneering Park C, Ningbo,315033,China.
Tel: +86-574-27709925 Fax:+86-574-27709939
Website:www.youjoy.cn E-mail:sales@youjoy.cn

EASE-01 Series Dental Unit

The ideal dental chair for what you want

Upgraded in 2023



Free, versatile, New-Fashion Choice

GIVE YOU MORE POSSIBILITIES

Ease 01 is a very high-quality dental chair among the same types on the market. It not only inherits the development concept of Runyes products, but also incorporates fashion elements. It is a trustworthy and ideal working partner for doctors in diagnosis and treatment.

The handpieces shown in the picture are options.

EASE-01U DENTAL UNIT

The multi-angle adjustable headrest can be adjusted quickly and easily with just one hand

The damping type hanging rack is clean and hygienic, greatly avoiding accidental equipment falling

Improved the structure of the air brake handle, which has a better grip and is easier to move the console.

Beautiful and comfortable sewing type leather

With linkage mechanism and chair lifting safety protection system

Covered Cushion Style Choices

- Standard domestic sofa leather(S01-S08)
- Optional American smooth seamless soft leather(I01-I13;N01-N02)
- Optional American Sewing Leather(I01-I13;N01-N02)



Streamlined and efficient operation process



Boutique and fashionable waist headrest (Optional, available for the same color with dental chair)



The side box cover can be opened on both sides. The water vapor and circuit system inside are isolated



Strong and weak suction filtration system



Disinfection system for options. Fulfill the needs of for disinfection of handpiece tubings, syringes, scaler.



Optional for surgical lamp



An assistant rack system that can be used to place objects or rotate at multiple angles



Lifetime Warranty Service for Water Pipe

Installation layout

