

EASE-02D DENTAL UNIT



Led light, hand control and induction dual control system provide more convenient operation

The unique luxury balance arm design highlights the product character and pursues the perfect combination of stability and smoothness

The micro switches are moved forward to the front, it is more convenient to operate

The hanger can be rotated 90 degrees to adjust the position according to the doctor's habits

Revision: 202302
Youjoy reserves the right to make modifications without prior notice.



Ningbo Youjoy Dental Equipment Co.,Ltd.
Add: No.275,Lane321,Tongning Road,Jiangbei Investment & Pioneering Park C, Ningbo,315033,China.
Tel: +86-574-27709925 Fax:+86-574-27709939
Website:www.youjoy.cn E-mail:sales@youjoy.cn



EASE-02 Series Dental Unit

Perfect fitting between appearance and function

Upgraded in 2023



Safe, Reliable, Comfortable

Ease-02's humanized safety design, strict production process, and full function make it easy for doctors and patients



Ergonomic design



Water purification system

New Ease 02 attracts many doctors with its elegant lines, unique shape, and powerful functions. For safety protection, the treatment instrument has a lifting safety protection device to ensure safe operation.

The handpiece shown in the picture is an optional accessory.

The multi-angle adjustable headrest can be adjusted quickly and easily with just one hand.

Covered Cushion Style Choices

- Standard American smooth leather
- Optional American sewing leather



With linkage mechanism and chair lifting safety protection system.

EASE-02U DENTAL UNIT



The damping type hanging rack is clean and hygienic, greatly avoiding accidental equipment falling.

Improved the structure of the air brake handle, which has a better grip and is easier to move the console.

Leather imported from American DCI.

Streamlined and efficient operation process



Water can be added directly without removing the water purification bottle.



Optional for surgical lamp.



An assistant rack system that can be used to place objects or rotate at multiple angles.



Disinfection system for options. Fulfill the needs of for disinfection of handpiece tubings, syringes, scaler.



PU base.



Lifetime Warranty Service for Water Pipe.



Boutique and fashionable waist headrest (Optional, available for the same color with dental chair).

Installation layout

