

Standard configuration: E96×2, E99×2, E70, E71×10

High configuration: E96×2, E99×2, E70×2, E71×10, E72×10

Compatible tips: E70, E71, E72, E73, E74, E93, E95, E96, E97, E98, E99

For canal irrigation after the root canal preparation

E70: Endochuck, used for install E71, E72, E73 or E74 for root canal cleaning or irrigation

E71, E72, E73, E74: Polymer material without cutting ability

E93: 20#, taper 1%, tip end length 21.5mm, Ni-Ti

E95: 20#, taper 2%, tip end length 18mm, Stainless steel

E96: 20#, taper 2%, tip end length 21mm, Stainless steel E97: 20#, taper 2%, tip end length 25mm, Stainless steel

E98: 25#, taper 2%, tip end length 18.5mm, Titanium alloy

E99: 25#, taper 2%, tip end length 18.5mm, Titanium alloy

Guilin Woodpecker Medical Instrument Co.,Ltd.
Information Industrial Park, Guilin National High-Tech
Zone, Guilin, Guangxi, 541004 P. R. China
Sales Dept.: +86-773-5873196/2350599
After-sales Service Dept.: 0773-5827898

E-mail: woodpecker4@glwoodpecker.com

Website: http://www.glwoodpecker.com zmn-xc-133 v1.0



WOODPECKER





