



Dental Turbine
HP33

LAST LONGER

SPECIFICATIONS

Model	HP33
Chuck	Push-button
Power	20W
Air pressure	2.5-2.7kgf/cm ²
Water pressure	0.5~2.0kgf/cm ²
Speed	≥280000r/min
Head size	Ø 12mm x H 13.1 mm

Guilin Woodpecker Medical Instrument Co.,Ltd.
 Information Industrial Park,Guilin National High-Tech
 Zone,Guilin,Guangxi,541004 P.R.China
 Sales Dept: +86-773-5873196/2350599
 After-sales Service Dept: 0773-5827898
 E-mail: woodpecker4@glwoodpecker.com
 Website: <http://www.glwoodpecker.com>

Scan and log on to the website
for more information



Find us on:



Find us on:




Find us on:



SERVICE LIFE

142.8%↑



Triple spray

Water-air separation
Effective cooling, reduced risk of burning tissues



German bearing

More stable




Titanium coating

Less wear

CLEAN HEAD SYSTEM

Prevent blood and foreign matter from entering the head,
enabling 100% zero suck-back
Avoid cross infection, extend service life

Well sealed clamping structure



Another brand

HP33: Zero suck-back

