


K-Po


Woodpecker brings you a comfortable treatment experience, farewell to the traditional polishing mode.

Wireless handpiece
Dual-control design
Cost-effective
Large-capacity lithium battery

Disposable prophy angle prevents cross-infection
Detachable and autoclavable outer sheath offers effective protection for motor handpiece



Dual control design enables easy operation
Wireless charging




Handpiece Control
Button control, Stable speed
6 gears optional:
500rpm 1000rpm 1500rpm
2000rpm 3000rpm 4000rpm

Foot pedal Control
Non-gear shift control mode
Speed range: 0 ~ 4000rpm
The speed will increase according to the depth of the pedal depressed


750mAh large capacity lithium battery
5h continuous working
A full charge can be used for at least 60 treatments

Noiseless high-quality motor

Small head brings greater comfort



Slim handpiece with anti-slip design



18mm diameter
Fish scale soft rubber design
guarantees non-slip grip.

① Prophy Angle Cup - Soft

Used for polishing tooth surface, preventing staining, perfecting and repairing tooth enamel

② Prophy Angle Cup - Hard

Used for polishing and cleaning teeth surface, bringing comfortable experience without harm to teeth, effectively removing biofilm and pigment

③ Prophy Angle Cone

Used for gingival nipple massage and adjacent surface polishing

⑤ Prophy Angle Brush - Sharp

Used to remove common stains on adjacent surfaces of teeth

④ Prophy Angle Brush - Flat

Used to remove common stains on the tooth surface

