

Surgical Loupes



Key Advantages




Selectable Angle of Lens

0° Lens



Simplified handling with direct vision.

48° Lens



Relieves cervical strain for lasting focus during extended procedures.



Adjustable Nose Support

Two options available for different nose bridge heights, ensuring a secure and comfortable fit.

Quick-Change Lens System


One frame accommodates multiple lenses, enabling easy switching without complex optical adjustments, significantly reducing usage and replacement costs.





Live Video and Projection Function

Built-in HD camera module for precise surgical documentation.



Supports real-time wireless projection for clinical teaching and team collaboration

Comfort-Fit Head Strap



- One-time adjustment for quick and easy daily use.
- Smooth buckle release prevents hair pulling and discomfort.

High-Definition Optics



Premium optical lenses deliver sharper details for more precise clinical outcomes.


Dental Headlight Series

The KP Dental Headlight Series offers both wired and wireless models, adaptable to different illumination levels and clinical environments, ensuring optimal visibility across diverse treatment scenarios.

i Light Wireless

Dental Headlight

Smart touch switch with multiple brightness levels, easy to operate.



Wireless freedom for greater flexibility and mobility in operation.




Equipped with two batteries, allowing quick replacement for continuous lighting.

i Light

Wired Headlight

Ultra-lightweight at just 24g, reducing neck and head fatigue for long procedures.




6-hour battery life ensures uninterrupted performance.



Compact lamp body with superior heat dissipation for stable, consistent illumination.

Surgical Loupes Specifications

Magnification	3.5X	4X	5X	3X	4X	5X
Weight	80g	84g	86g	82g	86g	89g
Angle of Declination	0° 			48° 		
Working Distance	320-420mm Loupe3532s	280-380mm Loupe4028s 360-460mm Loupe4036s 440-540mm Loupe4044s 500-600mm Loupe4050s	280-380mm Loupe5028s 360-460mm Loupe5036s 440-540mm Loupe5044s 500-600mm Loupe5050s	280-380mm Loupe3028 360-460mm Loupe3036 440-540mm Loupe3044 500-600mm Loupe3050	280-380mm Loupe4028 360-460mm Loupe4036 440-540mm Loupe4044 500-600mm Loupe4050	280-380mm Loupe5028 360-460mm Loupe5036 440-540mm Loupe5044 500-600mm Loupe5050

i Light Wireless Specifications

Headlight Unit	Weight	15g
	Light Diameter	Φ15mm
	Brightness	50000lux
	Color Temperature	5500K
	Power Consumption	3W
Battery	LED Lifespan	50000h
	Battery Capacity	300mAh
	Voltage	7.4V
	Runtime	1h
	Charging Time	1h
Charger	Input Voltage	110-240V
	Output Voltage	8.4V

i Light Specifications

Headlight Unit	Weight	24g
	Light Diameter	Φ18mm
	Brightness	70000lux
	Color Temperature	5500K
	Power Consumption	3W
Battery	LED Lifespan	50000h
	Battery Capacity	2800mAh
	Voltage	5V
	Runtime	6h
	Charging Time	4h
Charger	Input Voltage	100-240V
	Output Voltage	5V



**EQUIPOS E INSUMOS
ODONTOLÓGICOS**



**SERVICIO
TÉCNICO**



**8626-4641
8335-8987**



WWW.MEDIFIXCR.COM



VENTAS@MEDIFIXCR.COM

