



CE

ANYVAC SERIES ANYVAC 30

- Water saving system has reduced water consumption by up to 1/3 compared to existing products.
- Creates powerful vacuum pressure with less water
- Cover 3-4 unit chairs
- Easy to manage filter
- Voltage Min./Max. : AC 190 / 230V
- Size & Weight : 304×250×350mm, 18kg [120×98.5×138", 40lbs]

SPECIFICATION

| | | (as AC 220V) | |
|-------------------|-------|-----------------|-------|
| Frequency | Hz | 50 | 60 |
| Motor Output | kw | 0.62 | 0.75 |
| Full Load | A | 6.7 | 6.5 |
| Max. Vacuum | MPa | 0.065 | 0.070 |
| Max. Capacity | ℓ/min | 305 | 340 |
| Water Consumption | | 1.2 (as 0.6MPa) | |

ANYVAC SERIES ANYVAC40 (2h.p)

- Low noise with silencer-mounting type
- Prevents air pollution with water and air separation forms
- Efficient use of space
- Two pumps-mounting type to facilitate line change in case of malfunctioning
- Size & Weight : 520×425×830mm, 45kg [205×167×327", 99lbs]



ANYVAC SERIES SILENCER (Muffler)

- Improves health care environment by lowering noise
- Prevents air pollution with water and air separation forms
- Simple installation and management
- Size & Weight : 110×170×480mm, 3kg [43×67×189", 6.6lbs]

Made in MDMED (KOREA)