



MT2

Woodpecker
Dental Electric Motor

Refined trimming
Faster tooth preparation
More powerful crown-splitting



Imported German bearing
Low vibration enables
more precise tooth preparation

Imported German bearing
Stable and precise cutting
Guarantee stable speed and torque output

Minimally invasive repairing
precise teeth preparation

Compared with turbine handpiece, it has larger clamping force and can grasp the bur tighter. Effectively improve the precision of trimming and tooth preparation.

LCD screen
3.5 inches LCD screen
presents all functions.



Preparation with turbine handpiece
Rough surface with pits

Preparation with electric handpiece
Smooth surface without scratch

- Speed ratio
- Mode
- Speed adjustment
- LED
- Forward/Reverse rotation
- Setting
- Save

Parameters

Motor speed range	2000-40000rpm
(With our 1:5 speed-increase contra-angle, the maximum speed can reach 200000 RPM.)	
Maximum torque	5N.cm
Bur offset	<1.3um
Compatible contra-angle	16:1, 1:1, 1:5 (standard)
Motor size	φ22×76.7mm
Motor weight	70g

Configuration

- Main unit x1
- Motor with LED x1
- 1:5 contra-angle with LED and German bearing x1
- Power adapter x1
- 2-year warranty for main unit
- 1 year warranty for motor and contra-angle



High performance brushless motor
with greater torque and lower noise

Lower noise and faster tooth preparation.
Low abrasion, low heat generation and longer service life.



Four-hole water spray

Concentrated water output, better cooling effect.
Stable output effectively reduces aerosol, ensuring safer treatment environment.



One button control

More convenient and easier operation

Guilin Woodpecker Medical Instrument Co.,Ltd.
Information Industrial Park,Guilin National High-Tech
Zone,Guilin,Guangxi,541004 P.R.China
Sales department: 0773-5873196/2350599
Fax: 0773-5822490
E-mail: woodpecker4@glwoodpecker.com
Website: http://www.glzmn.com, http://www.glwoodpecker.com

