

I. Prevention Solution

PT-E

Introducing the Future of Dental Care



Dual Heating Systems

Provides warm water for both ultrasonic and airstream functions, enhancing patient comfort.

New Airstream Nozzle Design

Even water distribution, low-noise operation, and enhanced patient experience.

Advanced Airstream Handpiece with LED

Built-in LED for improved visibility and precision.

Titanium Alloy Handpiece

Lightweight, durable, and anti-slip for stable and comfortable grip.

Six Ultrasonic Modes:

Tailored for diverse clinical needs:

- G Mode: Ultra scaling, soft scaling
- P Mode: Periodontal treatment, implant maintenance
- E Mode: Endodontic treatment, irrigation

Illuminated Powder Container:

LED lighting differentiates powder types for easier identification, streamlining clinical operations:

White light: sodium bicarbonate powder

Green light: erythritol or glycine powder

