

# FAST

Woodpecker Prophylaxis Powder



## **FAST (PT-S1) : 65µm / 40µm Sodium Bicarbonate Supragingival (Lemon Flavor)**

For rapid removal of stubborn supragingival stains, calculus and biofilm in adults.

### Features

- 01** 300g/bottle
- 02** Sodium bicarbonate provides anti-inflammatory and anti-carries benefits
- 03** 2.5 Mohs' hardness ensures efficient and quick cleaning
- 04** Hardness between enamel and dentin, allowing for easy removal of supragingival biofilm and stubborn stains without damaging teeth