

FAST

Woodpecker Prophylaxis Powder



GENTLE MINI (PT-S3) : 25µm Glycine Supragingival + Subgingival

With a particle size of just 25µm, this glycine powder is gentle on soft tissues within periodontal pockets, gums, palate, and tooth surfaces. Suitable for pediatric dental care, orthodontic patients, restoration maintenance, periodontal care, and implant maintenance.

Features

- 01 **120g/bottle**
- 02 Fine particles for minimally invasive treatment
- 03 Hardness lower than tooth dentin, suitable for safe subgingival scaling, periodontal care, and implant maintenance

GENTLE (PT-S3) : 65µm Glycine Supragingival

For the removal of supragingival biofilm and soft deposits in adults. Suitable for patients choosing painless scaling, recall patients, sensitive patients and orthodontic patients.

Features

- 01 **200g/bottle**
- 02 Amino acids extracted from natural grains
- 03 Slightly sweet glycine with a pleasant taste
- 04 2.0 Mohs' hardness allowing gentle treatment